

Product Datasheet

Human FBLN1 ELISA Kit orb406244

Description	Human Fibulin-1(FBLN1) ELISA kit
Reactivity	Human
Range	6.25 µg/ml-400 µg/ml
Sensitivity	1.56 µg/ml
Tested Applications	ELISA
Target	FBLN1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm