

Product Datasheet

Human secretory component ELISA Kit orb406254

Description	Human secretory component,SC ELISA Kit
Reactivity	Human
Range	0.312 µg/ml-20 µg/ml
Sensitivity	0.078 µg/ml
Tested Applications	ELISA
Target	secretory component
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, urine, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm