

Product Datasheet

Human CYR61 ELISA Kit orb406271

Description	Human cysteine-rich 61 (Cyr61/CCN1) ELISA Kit
Reactivity	Human
Range	78 pg/mL-5000 pg/mL
Sensitivity	0.195 ng/mL
Tested Applications	ELISA
Target	CYR61
Note	For research use only.
Sample Types	serum, plasma, saliva, tissue homogenates
Uniprot ID	O00622
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, saliva, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm