

Product Datasheet

Human voltage-gated calcium channel antibody(IgG) ELISA Kit orb406285

Description	Human voltage-gated calcium channel antibody(IgG) ELISA Kit
Reactivity	Human
Range	Qualitative
Tested Applications	ELISA
Target	voltage-gated calcium channel antibody(IgG)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm