

## Product Datasheet

### Human IFNAR2 ELISA Kit orb406314

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Interferon alpha/beta receptor 2(IFNAR2) ELISA kit   |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 0.156 ng/ml-10 ng/ml   |
| <b>Sensitivity</b>         | 0.039 ng/ml  |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | IFNAR2   |
| <b>Note</b>                | For research use only.   |
| <b>Expiration Date</b>     | Please enquire.  |
| <b>Application Notes</b>   | <ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates, cell culture supernates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul> |