

Product Datasheet

Human PAPPA2 ELISA Kit orb406367

| | |
|----------------------------|---|
| Description | Human Pappalysin-2(PAPPA2) ELISA kit |
| Reactivity | Human |
| Range | 2 ng/ml -60 ng/ml |
| Sensitivity | 0.5 ng/ml |
| Tested Applications | ELISA |
| Target | PAPPA2 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma, urine, tissue homogenates.• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |