

Product Datasheet

Human PRDX6 ELISA Kit orb406391

Description	Human Peroxiredoxin-6(PRDX6) ELISA kit
Reactivity	Human
Range	12.5 pg/ml-800 pg/ml
Sensitivity	3.12 pg/ml.
Tested Applications	ELISA
Target	PRDX6
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm