

Product Datasheet

Human PPARG ELISA Kit orb406395

Description	Human Peroxisome Proliferator-activated receptor gamma;,PPAR-gamma; ELISA Kit
Reactivity	Human
Range	78 pg/ml-5000 pg/ml
Sensitivity	19.5 pg/ml
Tested Applications	ELISA
Target	PPARG
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm