

Product Datasheet

Human LMNA ELISA Kit orb406413

Description	Human lamin A/C (LMNA) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Tested Applications	ELISA
Target	LMNA
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P02545
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm