

Product Datasheet

Human MMP14 ELISA Kit orb406462

Description	Human Matrix Metalloproteinase 14, MMP-14 ELISA Kit
Reactivity	Human
Range	3.12 ng/ml-200 ng/ml
Sensitivity	0.78ng/ml
Tested Applications	ELISA
Target	MMP14
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm