

Product Datasheet

Human KLK6 ELISA Kit orb406481

Description	Human Kallikrein 6, KLK 6 ELISA Kit
Reactivity	Human
Range	1.25 ng/ml-80 ng/ml
Sensitivity	0.31 ng/ml
Tested Applications	ELISA
Target	KLK6
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm