

Product Datasheet

Human KLK6 ELISA Kit orb406481

Description	Human Kallikrein 6, KLK 6 ELISA Kit
Reactivity	Human
Range	1.25 ng/mL-80 ng/mL
Sensitivity	0.31 ng/mL
Tested Applications	ELISA
Target	KLK6
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q92876
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm