

Product Datasheet

Human PTH 1-34 ELISA Kit orb406494

Description	Human parathyroid hormone 1-34(PTH 1-34) ELISA kit
Reactivity	Human
Range	Request Information
Sensitivity	Request Information
Tested Applications	ELISA
Target	PTH 1-34 Ab
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm• Assay Principle: qualitative• Measurement: Indirect Experiment Notes: The antibody recognition range is 7-33 AA