

## Product Datasheet

### Human ARG2 ELISA Kit orb406538

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Arginase 2(ARG2) ELISA Kit  |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 0.47 mU/mL-30 mU/mL   |
| <b>Sensitivity</b>         | 0.12 mU/mL  |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | ARG2  |
| <b>Note</b>                | For research use only.  |
| <b>Sample Types</b>        | serum, plasma, cell culture supernates, urine, cerebrospinal fluid (CSF)  |
| <b>Uniprot ID</b>          | <b>P78540</b>   |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | • Sample Type: serum, plasma, cell culture supernates, urine and cerebrospinal fluid (CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm |