

Product Datasheet

Human TRAP ELISA Kit orb406542

Description	Human Tartrate Resistant Acid Phosphatase(TRAP)ELISA Kit
Reactivity	Human
Range	0.47 ng/ml-30 ng/ml
Sensitivity	0.11 ng/ml
Tested Applications	ELISA
Target	TRAP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm