

Product Datasheet

Human CSF1 ELISA Kit orb406544

Description	Human Macrophage Colony-Stimulating Factor, M-CSF ELISA kit
Reactivity	Human
Range	47 pg/ml-3000 pg/ml
Sensitivity	11.7pg/ml
Tested Applications	ELISA
Target	CSF1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm