

Product Datasheet

Human basal ganglia (BG) auto ELISA Kit orb406600

Description	Human anti basal ganglia (BG) auto ELISA kit
Reactivity	Human
Range	Request Information
Sensitivity	Request Information
Tested Applications	ELISA
Target	BG auto-Ab
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-3 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm