

Product Datasheet

Human streptolysin O antibody (IgG) ELISA Kit orb406615

Description	Human anti-streptolysin O antibody (IgG) ELISA Kit
Reactivity	Human
Range	Qualitative
Tested Applications	ELISA
Target	streptolysin O antibody (IgG)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm