

Product Datasheet

Cytokeratin 18 (Phospho-Ser52) antibody orb156372

Description	Rabbit polyclonal antibody to Cytokeratin 18
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from mouse CK18 around the phosphorylation site of Ser52: SR(p-S)VW
Form/Appearance	Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	P05783, Q5BJY9, P05784
Dilution Range	ELISA: 1:500-1000, IHC-P: 1:100-500, ICC: 1:100-500, IF: 1:100-500
Expiration Date	12 months from date of receipt.