

Product Datasheet

Human Transglutaminase 6 antibody (IgG) ELISA Kit orb406687

Description	Human Transglutaminase 6 antibody(IgG) ELISAKit
Reactivity	Human
Range	Request Information
Sensitivity	Request Information
Tested Applications	ELISA
Target	TGM6 Ab (IgG)
Note	For research use only.
Sample Types	serum
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm