

Product Datasheet

Human sCD146 ELISA Kit orb406698

Description	Human soluble cluster of differentiation 146,sCD146 ELISA Kit
Reactivity	Human
Range	1.25 ng/ml-80 ng/ml
Sensitivity	0.31 ng/ml.
Tested Applications	ELISA
Target	sCD146
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm