

Product Datasheet

Human CD28 ELISA Kit orb406702

Description	human soluble cluster of differentiation 28,sCD28 ELISA Kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml.
Tested Applications	ELISA
Target	CD28
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm