

Product Datasheet

Human IL7R ELISA Kit orb406711

Description	Human Soluble Interleukin-7 Receptor,sIL-7R ELISA Kit
Reactivity	Human
Range	0.78 ng/ml -50 ng/ml
Sensitivity	0.195 ng/ml
Tested Applications	ELISA
Target	IL7R
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm