

Product Datasheet

Human Galectin 3 ELISA Kit orb406712

Description	human lectin-galactose binding-soluble 3,Lgals3 ELISA Kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml
Tested Applications	ELISA
Target	Galectin 3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates, saliva• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm