

Product Datasheet

Human TK1 ELISA Kit orb406721

Description	Human TK1 ELISA Kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Tested Applications	ELISA
Target	TK1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P04183
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm