

## Product Datasheet

### Human sPECAM1/sCD31 ELISA Kit orb406724

<b>Description</b>	Human soluble platelet endothelial cell adhesion molecule 1,sPECAM-1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Sensitivity</b>	7.8 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	sPECAM1/sCD31
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Uniprot ID</b>	<b>P16284</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm