

Product Datasheet

Human TNFSF10 ELISA Kit orb406726

Description	Human soluble tumor necrosis factor-related apoptosis inducing ligand,sTRAIL ELISA KIT
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81 pg/ml
Tested Applications	ELISA
Target	TNFSF10
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, urine, cerebrospinal fluid(CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm