

Product Datasheet

Human HLA-B27 ELISA Kit orb406744

Description	Human leucocyte antigen B27,HLA-B27 ELISA Kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml
Tested Applications	ELISA
Target	HLA-B27
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, urine, cerebrospinal fluid (CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm