

Product Datasheet

Human HLA-G ELISA Kit orb406746

Description	Human leukocyte antigen G,HLA-G ELISA Kit
Reactivity	Human
Range	1.25 ng/ml-80 ng/ml
Sensitivity	0.31 ng/ml
Tested Applications	ELISA
Target	HLA-G
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm