

Product Datasheet

Human GM-CSF ELISA Kit orb406757

Description	Human anti-Granulocyte-Macrophage Colony Stimulating Factor(GM-CSF) antibody ELISA Kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml
Tested Applications	ELISA
Target	GM-CSF
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm