

Product Datasheet

Human PGAM2 ELISA Kit orb406766

Description	Human phosphoglycerate mutase 2,PGAM2 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9 pg/ml
Tested Applications	ELISA
Target	PGAM2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm