

Product Datasheet

Human PLA2R1 ELISA Kit orb406774

Description	Human phospholipase A2 receptor 1, 180kDa (PLA2R1) ELISA kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078ng/ml
Tested Applications	ELISA
Target	PLA2R1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm