

Product Datasheet

Human SERPINB3 ELISA Kit orb406785

Description	Human squamous cell carcinoma antigen 1 (SCCAg1) ELISA kit
Reactivity	Human
Range	0.31 ng/ml-10 ng/ml
Sensitivity	0.156 ng/ml
Tested Applications	ELISA
Target	SERPINB3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm