

Product Datasheet

Human HSPA5 ELISA Kit orb406831

| | |
|----------------------------|---|
| Description | Human immunoglobulin binding protein, BIP ELISA KIT |
| Reactivity | Human |
| Range | 0.156 ng/ml-10 ng/ml |
| Sensitivity | 0.039 ng/ml |
| Tested Applications | ELISA |
| Target | HSPA5 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |