

Product Datasheet

Human PAEP ELISA Kit orb406832

Description	Human Glycodelin ELISA Kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Sensitivity	0.195 ng/mL
Tested Applications	ELISA
Target	PAEP
Note	For research use only.
Sample Types	serum, cerebrospinal fluid (CSF), urine, tissue homogenates
Uniprot ID	P09466
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, cerebrospinal fluid (CSF), urine and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm