

Product Datasheet

Human PGRMC1 ELISA Kit orb406849

Description	Human Membrane-associated progesterone receptor component 1(PGRMC1) ELISA kit
Reactivity	Human
Range	47 pg/mL-3000 pg/mL
Sensitivity	11.7 pg/mL
Tested Applications	ELISA
Target	PGRMC1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O00264
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm