

Product Datasheet

Nm23H1 antibody {CF594} orb11140-CF594

Description	Rabbit polyclonal antibody against Nm23H1 conjugated to CF594.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine
Conjugation	CF594
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Nm23-H1: 41-152/152
Target	Nm23H1 {CF594}
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	17kDa
Uniprot ID	P15532, Q05982, P15531
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.