

Product Datasheet

Human LGALS3BP ELISA Kit orb406996

Description	Human Galectin-3-binding protein(LGALS3BP) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Sensitivity	0.39 ng/mL
Tested Applications	ELISA
Target	LGALS3BP
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q08380
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm