

Product Datasheet

Human ST antibody (IgM) ELISA Kit orb407022

Description	Human salmonella typhi (ST) antibody (IgM) ELISA kit
Reactivity	Human
Range	Qualitative
Tested Applications	ELISA
Target	ST antibody (IgM)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm