

Product Datasheet

Human C10orf59 ELISA Kit orb407059

| | |
|----------------------------|--|
| Description | Human Renalase/ Monoamine Oxidase-C(RNLS/ MAO-C)ELISA Kit |
| Reactivity | Human |
| Range | 0.625 µg/mL-40 µg/mL |
| Sensitivity | 0.156 µg/mL |
| Tested Applications | ELISA |
| Target | C10orf59 |
| Note | For research use only. |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | Q5VYX0 |
| Expiration Date | Please enquire. |
| Application Notes | • Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm |