

Product Datasheet

Human RBP4 ELISA Kit orb407079

Description	Human Retinol binding protein 4,RBP-4 ELISA Kit
Reactivity	Human
Range	15.6 ng/mL-1000 ng/mL
Sensitivity	3.9 ng/mL
Tested Applications	ELISA
Target	RBP4
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P02753
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm