

Product Datasheet

Human ALPP ELISA Kit orb407137

Description	Human Placental alkaline phosphatase,PLAP ELISA Kit
Reactivity	Human
Range	84 pg/ml-2000 pg/ml
Sensitivity	34 pg/ml
Tested Applications	ELISA
Target	ALPP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm