

Product Datasheet

Human GHbA1c ELISA Kit orb407155

Description	Human Glycated hemoglobin A1c, GHbA1c ELISA Kit
Reactivity	Human
Range	0.012 ng/mL-50 ng/mL
Sensitivity	0.012 ng/mL
Tested Applications	ELISA
Target	GHbA1c
Note	For research use only.
Sample Types	lysates for RBC
Expiration Date	Please enquire.
Application Notes	• Sample Type: Lysates for RBC • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm