

## Product Datasheet

# SARS-CoV-2 (2019-nCoV) Nucleocapsid (D63G, R203M, D377Y) Protein orb1173807

<b>Description</b>	Recombinant SARS-CoV-2 Nucleocapsid (D63G, R203M, D377Y) protein with C-terminal 6×His tag
<b>Target</b>	SARS
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Note</b>	For research use only.
<b>Purity</b>	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
<b>MW</b>	46.6 kDa
<b>Uniprot ID</b>	<b>P0DTC9</b>
<b>Expiration Date</b>	6 months from date of receipt.