

Product Datasheet

Human CNN1 ELISA Kit orb407169

Description	Human Calponin-1(CNN1)ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	CNN1
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	P51911
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm