

Product Datasheet

Human TSLP ELISA Kit orb407312

Description	Human thymic stromal lymphopoietin,TSLP ELISA Kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml
Tested Applications	ELISA
Target	TSLP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm