

Product Datasheet

Human PDGFC ELISA Kit orb407335

Description	Human vascular endothelial cell growth factor E,VEGF-E ELISA Kit
Reactivity	Human
Range	6.25 pg/ml-400 pg/ml
Sensitivity	1.56pg/ml
Tested Applications	ELISA
Target	PDGFC
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, urine, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm