

## Product Datasheet

### Human ANGPT2 ELISA Kit orb407338

<b>Description</b>	Human Angiopoietin 2,ANG-2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	39 pg/mL-2500 pg/mL
<b>Sensitivity</b>	9.75 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	ANGPT2
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>O15123</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm