

## Product Datasheet

### Human OxLDL ELISA Kit orb407396

<b>Description</b>	Human oxidized low density lipoprotein,OxLDL ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56 mU/mL-100 mU/mL
<b>Sensitivity</b>	0.78 mU/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	OxLDL
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, urine, saliva, tissue homogenates, cell lysates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• Sample Type: serum, plasma, urine, saliva, tissue homogenates, cell lysates.</li><li>• Assay Time: 1-5 hours</li><li>• Sample Volume: 50-100 <math>\mu</math>l</li><li>• Detection Wavelength: 450 nm</li></ul>