

Product Datasheet

Human IGF2BP2 ELISA Kit orb407411

Description	Human Insulin-like growth factor 2 mRNA-binding protein 2(IGF2BP2) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Tested Applications	ELISA
Target	IGF2BP2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y6M1
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm